



Be Right™



Beaker for Electrode/Sensor Calibration, colourless, 30 mL, pk/80

Številka izdelka: SM5010

EUR Cena: Kontakt

Na voljo

Colourless graduated 30 mL plastic beakers. Economic as it ensures optimum amount of buffer is used during calibration steps preventing from excess use of buffers by operators. Safe as disposable preventing from contamination. Can be used for both pH electrode and conductivity cell calibrations.

Obvezni pripomočki

- Beaker Holder (Item 923-556)
- Beaker Dispenser and Holder (Item 923-656)