



**Be Right™**



## **Beaker for pH Electrode Calibration, Green, 30 mL, pk/80**

**Številka izdelka:** SM5014

**EUR Cena:** Kontakt

**Na voljo**

Coloured and graduated 30 mL plastic beakers. It matches buffer color coding to prevent operators from wrong buffer use while calibrating. Economic as it ensures optimum amount of buffer is used during calibration steps preventing from excess use of buffers by operators. Safe as disposable preventing from contamination.

---

### **Obvezni pripomočki**

- Beaker Holder (Item 923-556)
- Beaker Dispenser and Holder (Item 923-656)