



Sodium Nitrate, 0.5M, 500 mL

Številka izdelka: 2507149

EUR Cena: Kontakt

Na voljo

Sodium Nitrate, 0.5 M, 500 mL, for 9245 and 9240 Sodium analysers

Specifikacije

Inštrument:	9245 & 9240
Koncentracija:	5 mg/L
Parameter:	Sodium
Pogoji skladiščenja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Prostornina / velikost paketa:	500 mL