



Kit of canisters for 9245-9240 sodium analyser.

Številka izdelka:	09240=A=8020
EUR Cena:	Kontakt
Odprema v 1 tednu	

Includes:

2 x bottles 0.5L all equipped with caps, stickers and tubing for Reactivtion and AutoCal

1 x Flask PolyEthylene 1L with anti-drop spout

Specifikacije

Inštrument:	9245 & 9240
Merilno območje:	N/A
Parameter:	Sodium
Pogoji skladiščenja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Prostornina / velikost paketa:	N/A