



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Številka izdelka: 6788302

EUR Cena: Kontakt

Na voljo

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Specifikacije

Komplet izdelka: 5500 sc Silica Analyser, 2/4 Channels

Vsebina paketa: Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Vsebina paketa

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions