



Chlorine Standard, 100-130 mg / L 20 mL Amp

Številka izdelka:	2835801
EUR Cena:	Kontakt
Na voljo	

This chlorine standard can be diluted down to verify the accuracy of the CL17 chlorine analyser or create a separate calibration curve. Under normal circumstances, a separate calibration curve is not required for the CL17. Package of one 20 mL ampule.

Specifikacije

Inštrument:	CL17
Merilno območje:	N/A
Pogoji skladiščenja:	Keep refrigerated.

Obvezni pripomočki

- Zaščitne rokavice velikosti L, modre, nitril, brez pudra, 50 kosov (Item SM743L)