



Konice pipet 1,0–5,0 mL

Številka izdelka:	LYW787
EUR Cena:	Kontakt
Na voljo	

Konice pipet 1,0–5,0 mL za pipeto s spremenljivo prostornino 1,0–5,0 mL, pipeto s stalno prostornino 0,2 mL in elektronsko pipeto 0,2–5,0 mL

Specifikacije

Opis:	Konice pipet (1,0–5,0 mL) za BBP065, BBP164 in BBP087
Parameter:	Konice za pipeto
Prostornina / velikost paketa:	1000 kosov