



## Kolorimeter POCKET II Komplet kolorimetra za kivetne teste za formaldehid

Številka izdelka: 5953000V.08

**ZASTAREL IZDELEK**  
Ta izdelek ni več na voljo.

Portable Colorimeter programmed for the determination of Formaldehyde with cuvette tests. In carrying case, with batteries and operating instructions.

Following Hach cuvette test can be evaluated with this instrument:

LCK325, 0.5-10.0 mg/L H<sub>2</sub>CO

Cuvette test has to be ordered separately.

**Simple: all functions are available via only four keys**

**Power: battery operation for a maximum of 2000 tests**

**Clear readings: even under difficult conditions, thanks to lightning and large figures on the display**

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifikacije

|                  |                                   |
|------------------|-----------------------------------|
| Merilno območje: | 0.5 - 10.0 mg/L H <sub>2</sub> CO |
| Parameter:       | Formaldehid (primerno za LCK325)  |