



Glass sample cell for lab turbidimeters, 30 mL, 6 pieces with caps

Številka izdelka:	2084900
EUR Cena:	Kontakt
Na voljo	

For use with 2100N, 2100AN, and Lab RATIO turbidimeters with caps

Specifikacije

Komentar:	For use with Hach Ratio 18900, 2100N and 2100AN Laboratory Turbidimeters.
Material:	Glass
Obseg dobave:	6 vials with caps
Prostornina vzorca:	30 mL
Teža:	0,2 kg