



HQ30D LBOD kit (cont. HQ30D, LBOD, probe, 6x BOD bottles, meter stand)

Številka izdelka:

HQ30D.99.701000

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Rugged portable meters for use in the field and plant.

The digital meter/electrode system combines reliability, flexibility and ease of use. Interchangeable Intellical electrodes are recognised automatically and store all relevant data. Rugged, nearly indestructible outdoor versions with multiple cable lengths enable measurements even at previously inaccessible locations.

Replaced by NEW HQ Series.

Specifikacije

Delovni pogoji: relativna vlažnost zraka:	90 % (non condensing)
Delovni pogoji: temperatura:	0 - 60 °C (32 - 140 °F)
Garancija:	36 mesecev
IP ocena ohišja:	IP67
Izvoz podatkov:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Merilno območje:	-1500 - 1500 mV
Meritev DO:	0.00 - 20.0 mg/L
Meritev slanosti:	0 - 42 g/kg
Meritev temperature:	-10 - 110 °C
Meritve prevodnosti:	0.01 µS/cm to 200 mS/cm
Merjenje upornosti:	2,5 Ωcm - 49 MΩcm
Model:	HQ30D – Multi/1 Channel
mV meritev :	-1500 to 1500 mV
mV resolucija:	0.1 mV
Natančnost prevodnosti:	± 0.5 in (1µS/cm - 200mS/cm) range
Needed Cable Length:	1
Parameter:	pH, mV, ISE, DO, Conductivity, TDS, Salinity, Resistivity, ORP, Temperature
Parameters:	BOD
	NA

	NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
pH Accuracy:	0,002 pH
pH Electrode calibration:	4 point calibration
pH Measurement Range:	0 - 14 pH
pH resolucija:	0.001 - 0.1 pH
Prikaz umeritvenih krivulj:	Pregled umerjanja
Resolucija:	0,5
Resolucija DO:	0.01 mg/L ali 0.1 %
Resolucija prevodnosti:	5 digits with 2 digits after decimal point
Sensor A:	LBOD10101 (1 m)
Sensor B:	NA
Sensor C:	NA
Sensors:	Sensor A: LBOD10101
	Sensor B: NA
	Sensor C:
Sestava:	No
Set?:	Da
Sonde vključene:	None
Spominska zmogljivost internih izmerjenih vrednosti:	500 results
TDS meritve:	0,0 - 50,0 mg/L
Temperaturna natančnost:	0.3 °C
Temperaturna resolucija:	0,1 °C
Teža:	0.323 kg
Type of measurement:	QC Lab
Umeritev DO elektrode:	0 - 1 DO calibration
Umeritveni intervali/alarmi/opomniki:	2 uri do 7 dni
Vhodi:	1
Vrsta sonde:	Intellical Digital SMART Standard Laboratory or Rugged Field
Vrsta zaslona:	240 ali 160 pixel
Vsebina paketa:	Meter package includes HQ30D Portable Meter; 4 AA batteries; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/ pk.
	Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.

Vsebina paketa

Meter package includes HQ30D Portable Meter; 4 AA batteries; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/pk. Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.