



Sension+ 9655 Fluoride Ion Selective Electrode (ISE) – sensor only

Številka izdelka:

LZW9655.97.0002

EUR Cena:

Kontakt

Na voljo

Sension+ 9655 is a half-cell (reference not integrated) fluoride ion selective electrode (ISE) with an epoxy body and non-replaceable crystalline membrane selective to fluoride ions in solution. It is recommended for use with the Sension+ 5044 reference electrode. The 9655 has a fixed 1 meter cable and BNC connector for laboratory use with the Hach Sension+ MM374 Multimeter. The 9655 is ideal for measuring fluoride concentrations in drinking water, wastewater and general water quality applications.

Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life or membrane replacement.

Provides fast, stable, and accurate response in a variety of sample types

Specifikacije

Dimenzija (d x g):	12 mm x 120 mm
Dolžina:	120 mm
Dolžina kabla:	1 m
Garancija:	6 mesecev
ISA pogoj:	258999
Material:	Sensor Body: Epoxy
Merilno območje:	0.02 mg/L - 20,000 mg/L F ⁻
Parameter:	Fluoride
Posebna značilnost:	Dry storage & fast response time. No replacement membranes.
Priključek:	BNC Connector
Raztopina za polnjenje:	-
Spoj sonde:	-
Temperaturno območje:	5 - 50 °C
Termistor:	-
Vrsta:	Lab ISE, potrebna referenčna
Vrsta senzorja:	Solid-state crystal membrane
Vrsta sonde:	Laboratory Half-Cell ISE
Vsebina paketa:	Includes: Sension+ 9655 Fluoride ISE (half-cell) and Basic User Manual.

Vsebina paketa

Includes: Sension+ 9655 Fluoride ISE (half-cell) and Basic User Manual.

Obvezni pripomočki

- Fluoride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 258999)
- Sension+ Electrolyte Solution 0.1M for ISE Electrodes, 125 mL (Item LZW9901.99)