



Sension+ ISE electrode, Nitrate (NO₃⁻), fixed cable 1m, BNC

Številka izdelka: LZW9662.97.0002

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Ion-selective indicator electrode for Nitrate measurements. Range 0.4-62,000 mg/L NO₃⁻, 5-40 °C. With 1m cable and BNC connector.

Specifikacije

Dimenzija (d x g):	12 mm x 120 mm
Dolžina:	120 mm
Dolžina kabla:	1 m
Garancija:	6 mesecev
ISA pogoj:	2984799
Material:	Sensor Body: Epoxy
Merilno območje:	0.1 mg/L (7x10 ⁻⁶ M) - 14,000 mg/L (1 M) NO ₃ ⁻ -N
Parameter:	Nitrate
Posebna značilnost:	Dry storage & fast response time. No replacement membranes.
Priključek:	BNC Connector
Raztopina za polnjenje:	-
Spoj sonde:	-
Temperaturno območje:	0 - 50 °C
Termistor:	-
Vrsta:	Lab ISE, potrebna referenčna
Vrsta senzorja:	Solid-state PVC membrane
Vrsta sonde:	Laboratory Half-Cell ISE

Obvezni pripomočki

- Nitrate Ionic Strength Adjustor Powder Pillows - Intellical, pk/100 (Item 2984799)
- Sension+ Electrolyte Solution 0.1M for ISE Electrodes, 125 mL (Item LZW9901.99)