



SENSION+ DO6 Field Kit with 5130 Polarographic DO Probe and case

Številka izdelka:

LPV4552.98.0002

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Complete mobile measurement stations.

Kit contains DO portable meter, polarographic DO probe, portable field carrying case and accessories.

The basic Dissolved Oxygen Field Kit is an all-in-one portable laboratory that makes general water quality meter testing in the field fast and simple. The field kit is designed to be used in a wide variety of applications and comes complete with everything needed to start testing.

The kit is ideal for dissolved oxygen measurement in general aqueous applications – particularly in a field environment.

Simple and fast measurements

Ergonomic design for real one hand operation

IP67 protected

Precise and cost-efficient system.

Specifikacije

AC in USB delovanje:	No
Delovni pogoji: relativna vlažnost zraka:	< 80 % relative humidity (non-condensing)
Delovni pogoji: temperatura:	0 - 50 °C (32 - 122 °F)
Delovni vmesnik:	Keypad
Dimenzija (v x š x g):	38 mm x 74 mm x 186 mm
Garancija:	36 mesecev
GLP značilnosti:	-
IP ocena ohišja:	IP67
Izhod:	-
jeziki:	-
Jeziki v navodilih za uporabo:	Multi-lingual
Merilna metoda:	Auto stabilisation, continuous
Merilno območje:	-
Meritev barometričnega tlaka:	Manual introduction
Meritev DO:	0.00 - 19.99 mg/L and 20.0 - 22.0 mg/L (25°C)

Meritev slanosti:	-
Meritev temperature:	-20.0 - 150.0 °C
	(-4.0 - 302.0 °F)
Meritve prevodnosti:	-
Meritve prevodnosti pri stabilnem branju:	-
Meritve prevodnosti: popravek temperature:	-
Model:	DO6
mV resolucija:	-
Namizni:	No
Natančnost:	-
Natančnost prevodnosti:	-
Osvetlitev ozadja:	Yes
Parameter:	Dissolved Oxygen, Temperature
PC data transfer software :	-
pH Accuracy:	-
pH Electrode calibration:	-
pH Measurement Range:	-
pH resolucija:	-
Potrdila o skladnosti:	CE.WEEE
Pressure compensation:	Manual input
Resolucija DO:	0.01 mg/L depending on measurement range
Resolucija prevodnosti:	-
Resolution:	0.01 mg/L
Salinity resolution:	-
Samodejno prepoznavanje pufra:	-
Set?:	Da
Simultaneous measurements:	Dissolved Oxygen, Temperature
Sonde vključene:	5130 Electrode
Spominska zmogljivost internih izmerjenih vrednosti:	-
TDS meritve:	-
TDS resolucija:	-
Temperaturna kompenzacija:	Automatic or manual
Temperaturna natančnost:	0,2 °C (± 1 številka)
Temperaturna resolucija:	0,1 °C
Teža:	0.3 kg
Umeritev DO elektrode:	1 to 2 point, saturated air, 0% DO calibration chamber delivered with kit
Upor vode:	IP67
Vhodi:	MP-5
Višina (skupaj):	-
Vmesnik: jeziki:	Icon based, no language needed
Vrsta zaslona:	128 x 64 pixel, simultaneous readings
Zahteve baterije:	3, AA
Zaslon: LCD (osvetlitev iz ozadja):	Backlit LCD with pictograms

Obvezni pripomočki

- Filling Solution, Internal, Polarographic DO, 30 mL (Item 2759123)
- Filling Solution, Internal, Polarographic DO, 50 mL (Item 2759126)