



PHC3359 Combined pH electrode for micro samples, fibre, BNC

Številka izdelka: **E16M343**
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

pH electrode with Ag/AgCl reference system. Measuring range 0-12 pH, 0-80°C. Liquid junction: fibre. With BNC plug.

Specifikacije

Garancija:	6 mesecev
Parameter:	pH
Priključek:	BNC Connector
Vrsta:	Lab, kombinirana elektroda
Združljivi instrumenti:	MeterLab & sensION+ GLP lab