



Titration Application Pack Water: pH, Alkalinity & Conductivity

Številka izdelka: **AP0002.AT1102**

EUR Cena:
Odprema v 2 tednih

Kontakt

Application Pack for the determination of pH, Alkalinity and Conductivity in water. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Aplikacija:	Municipal Applications
Garancija:	24 mesecev
Parameter:	pH/Alkalinity & Conductivity
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Združljivi instrumenti:	AT1102
	AT1112
	AT1122
	AT1222

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.