



Titration Application Pack Karl Fischer: Moisture Content

Številka izdelka: **AP0014.KF1121**

EUR Cena: Kontakt

Odprema v 2 tednih

Application Pack for the determination of Moisture content (Karl Fischer). Suited for Titrallab KF1000 Series, Automatic Titrator Model KF1121.

Specifikacije

Aplikacija:	Moisture Content Application
Garancija:	24 mesecev
Parameter:	Moisture Content (Karl Fischer)
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Združljivi instrumenti:	KF1121

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.