



TitraLab KF1000 Series Vol. Karl Fischer Titrator, 1 Syringe, 2 Pumps

Številka izdelka: **KF1121.98**

EUR Cena: Kontakt

Odprema v 3 tednih

Simple and accurate volumetric titrations to measure moisture content

Simplify complex volumetric Karl Fischer titration with the Titralab KF1000 from Hach. Straight from the box the KF1000 to easily and quickly perform moisture content determinations. The automatic system gives you more reliable results without complicated analysis. To make setup and titration simple, the KF1000 utilizes the application specific kit to make it quick and easy for anyone to set up and operate a test.

Specifikacije

Aplikacija:	5
Bireta: možni volumni:	2.5/5/10/25 mL, ISO8655-3
Burette for Titrant:	For Titrant 1
Certificiranje:	Safety IEC/EN 61010-1; EMC IEC/EN 61326-1
Črpalke:	2
Delovna temperatura:	15 - 35 °C
Delovna vlažnost:	20 - 80 %
Elektroda:	MTC Intellical ("Plug & Play" Digital)
Frekvenca (Hz):	50/60 Hz
Garancija:	24 mesecev
Izvoz podatkov:	On USB flash drive, CSV format, Excel compatible
Material:	PP Latene, silicone, hardened glass, POM, stainless steel. Splashproof design.
Merilna metoda:	Volumetric Karl Fischer, AC or DC imposed current
Merilni način:	Sample, blank, sample with blank, QC sample, QC sample with blank
Model:	KF1121
Napetost (V):	100/240 V AC
Pogoji skladiščenja:	-5 - 40 °C
Pomnilnik/hranjenje podatkov:	Last 100 samples, QC & blank analyses, last 10 calibrations
Prikaz umeritvenih krivulj:	KF Titration curve
Pump flow rate:	170 mL/min
Required Application Package/s:	AP0014.KF1121 Moisture Content

Samodejno dodajanje reagenta:	2 integrated peristaltic pumps for reagent addition
Stojalo za vzorec:	Integrated, magnetic stirring, KF glass beakers 150 mL
Supported peripherals:	Printer, PC software, balance, keyboard, mouse, barcode reader
Teža:	4 kg
Umeritev:	Titrant
Umeritveni intervali/alarmi/opomniki:	Da, uporabniško določeno obdobje
Volumen čaše sprejet:	System specific beakers: Proprietary glass measurement cell, screw-on
Vrsta:	Karl Fischer
Vsebina paketa:	1 ea. KF1121-Titralab KF1000 Titrator
	1 ea. LZE144- Power Supply
	1 ea. LZE141- 4 position tube holder
	1 ea. Basic User Manual
	Sold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.
Zaščita z geslom:	Yes, user defined

Vsebina paketa

1 ea. KF1121-Titralab KF1000 Titrator 1 ea. LZE144- Power Supply 1 ea. LZE141- 4 position tube holder 1 ea. Basic User ManualSold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.

Obvezni pripomočki

- Titration Application Pack Karl Fischer: Moisture Content (Item AP0014.KF1121)