



Nitrate standard, 200 mg/L NO₃ (45.2 mg/L NO₃-N), 500 mL

Številka izdelka: LCW827

EUR Cena: Kontakt

Na voljo

Standard solution for online nitrate determination

Specifikacije

Inštrument:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Koncentracija:	200 mg/L NO ₃ (45.2 mg/L NO ₃ -N)
Parameter:	Nitrate
Platforma:	Raztopina
Prostornina / velikost paketa:	500 mL