



## EDTA Standard Solution, 0.200 N, 500 mL

Številka izdelka:

102149

Nearno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



---

### Specifikacije

|                      |                                                                                                 |
|----------------------|-------------------------------------------------------------------------------------------------|
| Inštrument:          | SP510                                                                                           |
| Koncentracija:       | 0.200 N                                                                                         |
| Merilno območje:     | 1000 - 25,000 mg/L CaCO <sub>3</sub>                                                            |
| Parameter:           | Hardness                                                                                        |
| Platforma:           | Raztopina                                                                                       |
| Pogoji skladiščenja: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |
| Število testov:      | 40 - 100                                                                                        |