



Beaker, polypropylene, low form, 600 mL

Številka izdelka: 108052

EUR Cena: Kontakt

Odprema v 2 tednih

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Specifikacije

Enota:	1 each
Inštrument:	All except DR/4000, PC II
Kapaciteta:	600 mL
Material:	Polypropylene, Low Form
Premer, zunanji:	110.7 mm with spout
Številka metode:	8006