



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Številka izdelka: 108076

EUR Cena: Kontakt

Odprema v 2 tednih

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Enota:	6 /pk
Kapaciteta:	250 mL
Material:	Polypropylene, Low Form
Material ohišja:	Polypropylene
Premer, zunanji:	83.6 mm with spout