



**Be Right™**



## **Funnel, analytical, 114 mL approx. volume**

**Številka izdelka:** 108368

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### **Specifikacije**

Prostornina / velikost paketa: 1 each