



## Buffer Solution, pH 7.00, 25 mL

Številka izdelka: 1222220

EUR Cena: Kontakt

Odprema v 2 tednih

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colourless pH buffer solution. pH is  $7.00 \pm 0.02$  at 25°C. 25 mL bottle

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

Barva kodirana?:	Ne
Opis:	pH Buffer Calibration reagent
Package settings:	Bottle
pH vrednost(i):	$7.00 \pm 0.02$ pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf-life after opening, following recommended storage and use: 6 months
Prostornina / velikost paketa:	25 mL
Sledljivo po NIST-u:	Da