



Flask, volumetric, polypropylene, 100 mL

Številka izdelka: 1406042

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifikacije

Certificiranje:	Not Applicable
Enota:	1 each
Kapaciteta:	100 mL
Material ohišja:	Polypropylene