



Flask, volumetric, polypropylene, 250 mL

Številka izdelka: 1406046

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifikacije

Certificiranje:	Not Applicable
Enota:	1 each
Kapaciteta:	250 mL
Material ohišja:	Polypropylene