



Test kit, nitrate-nitrite, model NI-12

Številka izdelka:

1408100

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nevarno



Simple, low-cost color disc method is the most accurate visual test kit method.

Specifikacije

Ime metode:	Color disc/Cadmium reduction, Diazotization
Merilno območje:	0 - 40 mg/L NO_3^- -N
Model:	NI-12
Parameter:	Nitrogen, Nitrate/Nitrite
Platforma:	Color Disc
Smallest increment:	1 0.01
Standardno merilno območje:	0 - 40 mg/L NO_3^- -N; 0.01 - 0.40 mg/L NO_2^- -N
Število testov:	100
Test Kit model:	NI-12**
Vrsta primera:	D