



Test kit, copper (free) test color disc, model CU-5

Številka izdelka:

1421300

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. This kit is USEPA-accepted for waste water and potable water and approved by some states for potable water under specific conditions (Check with your regulatory agency). Range: 0.1-5 mg/L; contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color-disc/Bicinchoninate
Merilno območje:	0 - 4 mg/L Cu
Model:	CU-5
Parameter:	Copper (free)
Platforma:	Color Disc
Smallest increment:	0.1
Standardno merilno območje:	0 - 4 mg/L Cu
Število testov:	100 F
Test Kit model:	CU-5
Vrsta primera:	B