



## Standard solution, chloride, 12,500 mg/L as Cl<sup>-</sup> (NIST), pk/16 - 10 mL voluette® ampules

Številka izdelka: 1425010

EUR Cena: Kontakt

Odprema v 2 tednih

For accuracy checks in chloride determinations

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Specifikacije

Koncentracija:	12500 mg/L Cl <sup>-</sup>
Opis:	Chloride Standard Solution as Cl <sup>-</sup> (NIST)
Platforma:	Raztopina