



Standard solution, iron, 25 mg/L as Fe (NIST), pk/16 - 10 mL voluette® ampules

Številka izdelka:

1425310

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For accuracy checks in iron determinations

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	25 mg/L
Opis:	Iron Standard Solution as Fe (NIST)
Platforma:	Raztopina
Prostornina / velikost paketa:	mL/16