



Volatile Acid Standard Solution, 62,500 mg/L as Acetic Acid, pk/16 - 10 mL Voluette Ampules

Številka izdelka: 1427010

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Acetate Standard Solution. For accuracy checks in volatile acids determinations.

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	62500 mg/L
Platforma:	Raztopina
Prostornina / velikost paketa:	10 mL/ pack of 16