



Manganese test kit, model MN-5

Številka izdelka: **146700**
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/Cold periodate
Merilno območje:	0 - 3.0 mg/L Mn
Model:	MN-5
Parameter:	Manganese
Platforma:	Color Disc
Smallest increment:	0.1
Standardno merilno območje:	0 - 3.0 mg/L Mn
Število testov:	100
Test Kit model:	MN-5
Vrsta primera:	D