



Acetate Buffer Solution, pH 4.0, 100 mL

Številka izdelka:

1490932

Nearno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



Acetate Buffer Solution, pH 4.0. Total and/or Free Chlorine. Method: Amperometric. pk/100 mL.

Specifikacije

Ime metode:	Amperometric Titration
Inštrument:	AutoCat [™] 9000
Merilno območje:	6 - 1000 µg/L Cl ₂
Platforma:	Raztopina
Številka metode:	8209
Število testov:	100