



Hardness, iron, and pH test kit, model HA-62B

Številka izdelka:

183702

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Model HA-62B

Includes color disc and color comparator for colorimetric comparison of iron and pH, as well as a drop-count titrator for hardness determination. Test kit produces a wide range of pH tests.

Combines frequently measured parameters in one rugged kit

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Aplikacija:	Water Quality
Model:	HA-62B
Parameter:	Hardness, Iron, and pH (4-10)
Platforma:	Color Disc
Število testov:	100