



Hydrazine color disc test kit, model HY-2

Številka izdelka:

184900

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/ p-Dimethylamino- benzaldehyde
Merilno območje:	0 - 1 mg/L N_2H_4
Model:	HY-2
Parameter:	Detergents
Platforma:	Color Disc
Smallest increment:	0
Standardno merilno območje:	0 - 1 mg/L N_2H_4
Število testov:	300
Test Kit model:	HY-2
Vrsta primera:	D