



Cyanuric acid test kit, model CY-3

Številka izdelka: 185102

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Aplikacija:	Pool and Spa
Ime metode:	Turbidimetry
Merilno območje:	20 - 100 mg/L
Model:	CY-3
Parameter:	Cyanuric Acid
Smallest increment:	20
Število testov:	100
Test Kit model:	CY-3
Vrsta primera:	D