



Cylinder, graduated, mixing, glass, 25 mL ± 0.3 mL, 0.5 mL divisions, glass stopper #13

Številka izdelka: 189640

EUR Cena: Kontakt

Odprema v 2 tednih

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	25 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Vrsta:	eno metričen, mešan