



## Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Številka izdelka: 190902

EUR Cena: Kontakt

Odprema v 2 tednih

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

---

### Specifikacije

|                                |                    |
|--------------------------------|--------------------|
| Body Material:                 | Borosilicate Glass |
| Kapaciteta:                    | 60 mL              |
| Material:                      | Borosilicate Glass |
| Package size UOM:              | 1pk                |
| Prostornina / velikost paketa: | 1 each             |