



Sodium hydroxide standard solution, 0.020 N (N/50), 1 L

Številka izdelka:
EUR Cena:
Odprema v 2 tednih

19353
Kontakt

Nearno



For acidity determination by buret titration. Hach methods 8219 (methyl orange) and 8010 (phenolphthalein). 1 liter bottle

Specifikacije

EPA compliant:	Da
Koncentracija:	0.020 N
Merilno območje:	0 - 10000 mg/L CaCO ₃ (P)
Obseg:	1 L
Platforma:	Raztopina
Številka metode:	8288 (Acid)