



Be Right™



Lactose broth, dehydrated 500 g

Številka izdelka: 19634
EUR Cena: Kontakt
Pokliči za datum odpreme

Hach microbiological chemicals

Hach carries a wide variety of media and accessories for microbiology applications, including dehydrated media, nutrient agar, broths, acids, dilution water, sterilization indicator, and reagents. Combined with Hach's microbiology expertise, these products mean making your own media in the lab for standard or custom sampling applications is economical and practical. Instructions for mixing are included on each bottle. Dehydrated media product is available for your application, including m-Endo, Nutrient, Tryptic, m-Green, and many others.

Specifikacije

Media Type:	Lactose broth
Prostornina / velikost paketa:	500 g
Region:	Global
Volume / package size e Post Text:	500 g