



Test kit, hardness & iron color disc, model HA-77

Številka izdelka: 202300
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

The kit to use for more resolution in a visual test

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Aplikacija:	Water Quality
Model:	HA-77
Parameter:	Hardness and Iron
Platforma:	Color Disc