



Test kit, hardness (total), model HA-DT

Številka izdelka: 2063600

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Utilizing the Hach Digital Titrator, the Model HA-DT kit is used to determine hardness concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikacije

Ime metode:	Digital Titrator/EDTA
Merilno območje:	10 - 4,000 CaCO ₃
Najmanjši korak rasti:	0.1-10
Parameter:	Hardness, total - As CaCO ₃
Platforma:	Digitalni Titrator
Število testov:	100
Test Kit model:	HA-DT
Vrsta primera:	F