



Test kit, carbon dioxide, model CA-DT

Številka izdelka:

2064100

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Utilizing the Hach Digital Titrator, the Model CA-DT kit is used to determine carbon dioxide concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10.1-1,000 mg/L. Contains approximately 50 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikacije

Ime metode:	Digital Titrator/ Sodium hydroxide Phenolphthalein
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Merilno območje:	10 - 1,000 mg/L CO ₂
Parameter:	Carbon Dioxide - As CO ₂
Platforma:	Digitalni Titrator
Smallest increment:	0.1, 2.0
Število testov:	50
Test Kit model:	CA-DT
Vrsta primera:	F