



Cone, Imhoff settling, 1 L, polycarbonate, w/ screw closure

Številka izdelka: 206700

EUR Cena: Kontakt

Odprema v 2 tednih

Fully graduated cones for determination of small amounts of settleable matter in sewage by Imhoff Sedimentation Method. Glass cone features reinforced tip for ruggedness, while the Polycarbonate cone features a leak-proof, removable closure for withdrawal of sediment for analysis. Support ring, support stand, and cone holders available and sold separately.

Specifikacije

Kapaciteta:	1 L
Prostornina / velikost paketa:	1 each