



Be Right™



Tubing, nalgene 8000, vacuum, 9.5 mm (I.D.)

Številka izdelka: 2074146

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Similar to Nalgene 8000 regular tubing, but with heavier walls. Imprinted in red for instant identification. Saves time and avoids errors which could invalidate experiments or injure laboratory personnel.

Specifikacije

Debelina stene cevi:	1/4 in
Dolžina:	50 ft
Formula (Text):	R-3603
I. D.:	3/8 in
Material ohišja:	Tygon
Premer, zunanji:	7/8 in