



Be Right™



Cylinder, graduated, mixing, Glass, 25 mL ± 0.3 mL, 0.5 mL divisions, glass stopper #13

Številka izdelka: 2088640

EUR Cena: Kontakt

Na voljo

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	25 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Vrsta:	eno metričen, mešan