



## Flask, Erlenmeyer, glass w/screw cap, 250 mL

Številka izdelka: 2089646

**ZASTAREL IZDELEK**  
Ta izdelek ni več na voljo.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### Specifikacije

|             |                 |
|-------------|-----------------|
| Enota:      | 1 each          |
| Kapaciteta: | 250 mL          |
| Material:   | Glass Screw Cap |