



Flask, Erlenmeyer, glass w/screw cap, 500 mL

Številka izdelka: 2089649

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Predlagani nadomestki:
• 2089679

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specifikacije

Enota:	1 each
Kapaciteta:	500 mL
Material:	Glass Screw Cap