



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Številka izdelka: 2089676

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specifikacije

Enota:	6/pk
Kapaciteta:	250 mL
Material:	Glass Screw Cap