



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 500 mL

Številka izdelka: 2089679

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specifikacije

Enota:	6/pk
Kapaciteta:	500 mL
Material:	Glass Screw Cap