



Be Right™



Flask, Erlenmeyer, Polycarbonate Capacity 125 mL

Številka izdelka: 2089743

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Graduated (Molded-in graduations) flask features: clarity, toughness and strength. Autoclavable; however, repeated autoclaving may result in loss of mechanical strength. #5 Neoprene stoppers (P/N 14808-05; pk/12) available separately. 125 mL.

Specifikacije

Enota:	1 each
Kapaciteta:	125 mL
Material:	Polycarbonate